

Nitrous oxide Analysis

Separation of CO₂ & N₂O (m/z 44) by using infiTOF

Introduction

Nitrous oxide (N₂O) is known as a greenhouse gas, and the warming effect is about 300 times larger than carbon dioxide (CO₂). Furthermore, N₂O is one of the ozone-depleting substances. So it is quite significant to monitor N₂O.

To carry out real-time monitoring of N₂O using mass spectrometry, a mass spectrometer with high mass resolution is required, because the nominal mass of N₂O is the same as that of CO₂. To separate CO₂ and N₂O doublet completely, the required mass resolution is about 10,000. Conventional quadrupole mass spectrometers (QMS) are incapable of this measurement because of low mass resolution.



infiTOF

The infiTOF system is compact, portable, and capable of achieving a resolving power of 30,000. High-resolution mass analysis of CO₂ and N₂O with infiTOF is possible.

Sample gas; Carbon dioxide (CO₂) and Nitrous oxide (N₂O) mixture gas

MS conditions; Ion source : EI(Pos), Ionization voltage : 70eV, Ion source Temperature : 250°C

Sample gas was introduced directly to the ion source through the fused silica capillary tube.




Fig. 1 Separation of CO₂ and N₂O doublet. (A) 2cycles, (B) 10cycles, (C) 20cycles and (D) 50cycles. At 100 cycles, mass resolution of 22,000 was achieved.

Conclusion

In the short flight length (2cycles, Fig.1A), the doublet was still unified. After 10 cycles, the doublet began to separate, but not yet completely separated (Fig. 1B). After 20 cycles, CO₂ and N₂O were clearly separated, and the obtained mass resolution was 11,000 (Fig.1C). After 50 cycles, the mass resolution of 22,000 was achieved (Fig.1D).



KANOMAX
 The Ultimate Measurements

KANOMAX CORPORATION

<Tokyo Office>

2-6-2 Hamamatsu-cho, Minato-ku, Tokyo 105-0013 JAPAN

TEL: +81-3-5733-6544

FAX: +81-3-5733-6545

E-mail: analytical@kanomax.co.jp

URL: <https://kanomax.biz/asia/index.html>

AAN002-2E