

# Microbial Air Sampler

Model 3080/3081



# User-Friendly Next Generation Microbial Air Sampler

### **Features**


- Simple operation with 4.3inch touch screen
- High precision automatic flow control
- Protecting data against unauthorized access by administrative settings
- SUS316 sampling head is adopted to prevent contamination by microbial
- USB / RS-485 communication, operation by remote controller
- External printer available






#### Specification **Product Name** Microbial Air Sampler Model Number 3080 3081 Sampling Method Inertial Impaction Method 100L/min (Accuracy: $\pm 2.5\%$ ) Flow Rate **Exhaust Filter** O(Class ISO3) X Sampling Volume 20~6.000 L Sampling Mode Single, Repeat (multi-cycle sampling) Sampling Times $2\sim300$ times, or continuous sampling $10 \text{sec} \sim 1 \text{ hour (when sampling mode is "repeat")}$ Sampling Interval Applicable Culture Dish Petri dish(outer diameter: 90~100mm)/Contact plate(outer diameter: 55~84mm) Material Sampling head case: SUS316, Front panel and handle: ABS resin Administrator / General user authorization settings, user name, password settings Administrative Settings 4.3inch touch screen, Sampling volume and sampling time can be displayed **Operation Function** Language can be switched either Japanese, English or Chinese Printer External printer Communication Function USB / RS-485 Memory In-built memory Storage Method Data storage and Method Max. 10,000 data, CSV format Internal power lithium-ion battery AC Adopter (AC100~240V) Outer power Power Continuous battery operation time Max. 7 hour 2~3hours Battery recharging time **Operation Environment** Temperature: $0\sim40^{\circ}$ C, Humidity: $0\sim90\%$ (\*No condensation) **Dimensions** W136 $\times$ H235 $\times$ D133 mm Weight Approx. 3.2 kg 3.0 kg

### Dimension of unit







Carrying case Model 3080-70



## KANOMAX JAPAN, INC.

2-1, Shimizu, Suita-shi Osaka, 565-0805

TEL: 81-6-6877-0183 FAX: 81-6-6877-5570 E-mail: sales@kanomax.co.jp URL: www.kanomax.co.jp www.kanomax-usa.com Distributed by :